## Scanning Tunnelling Microscopy (S.T.M)

A Scanning Tunnelling Microscope is an instrument for imaging surfaces at the atomic level. It is development in 1981 earned its inventors, Gerd Binning and Heinrich Rohrer (at IBM Zürich), the Nobel Prize in Physics in 1986. For an STM, good resolvhon is considered to be 0.1 nm lateral resolvhon and 0.01 nm (10 pm) depth resolvhon. With this resolvhon, individual atoms within materials are routinely imaged and manipulated. The STM can be used not only in ultrating vacuum but also in air, water and various other liquid or gas ambients, and at temperatures ranging from OX to over 1000°C.

## Working Principle

Scanning Tunnelling Microscope's principle of operation is based on the quantum mechanical phenomenon known as hindling, in which the wavelike propertie of electrons permit them to "tunnel" beyond the surface of a solid into regions of space that are forbidden to them under the rules of classical physics. The probability Of finding such tunnelling electrons decreases exponentially as the distance from the surface increases. The STM makes use of this extremely sensitivity to distance. The sharp tip of a tungsten needle is positioned a few angstroms from the sample surface. A small voltage is applied between the probe hip and the surface, Causing the electrons to tunnel across the gap. As the probe is scanned over the surface, it registers Variations in the tunnelling current, as the and this information can be processed to provide a topographical image of the surface. Tunnelling is a functional functioning concept that

arises from quantum mechanics. Classically, an object hithing an impenetrable barrier will not pass through. In contrast, objects with a very small mass, such a, the electron, have wavelike characteristics which permit such an event, referred to as tunnelling.

Electrons behave as beams of energy, and in the presence of a potential U(z), assuming 1-dimensional case, the energy levels  $\Psi_n(z)$  of the electrons are given by Solutions to Schrödingers equation.

$$-\frac{\hbar^{2}}{2m}$$
.  $\frac{\partial^{2} \Psi_{n}(z)}{\partial z^{2}} + U(z) \Psi_{n}(z) = E \Psi_{n}(z)$ 

where to is the reduced Planck's constant, z is the position, and m is the mass of an electron. If an electron of energy E is incident upon an energy borrier of height U(z), the electron wave function is a travelling wave solution,

$$\psi_n(z) = \psi_n(0) e^{\pm ikz}$$
where  $k = \sqrt{2m(E - U(z))}$ 

If E > U(z), which is true for a wave function inside the hip or inside the sample Inside a barrier, inside the hip or inside the sample Inside a barrier, inside the hip or inside the sample Inside a barrier, inside the hip or inside the sample Inside a barrier, and  $E \neq U(z)$  so the wave functions which satisfy this are decaying waves,  $\forall n(z) = \forall n(o) \in \pm kz$ 

whore k= 
$$\sqrt{2m(U-E)}$$

quantities the decay of the wave inside the barrier, with the barrier in the tz direction of K.

1950384

1862, 9

## CONSTRUCTION OF STANDARD TUNIFELLING MICROSCOPE

First,

close

cor

The components of an STM include Scanning tip, piezo electric controlled height and X, x scanner coarse sample-to-tip control, Vibration isolation system, and a computer.

The resolution of an image is limited by the rodus of curvature of the Scanning hip of the STM.

Additionally, image artifacts can occur if the hip has 2 tips at the end rather than a single atom; this leads to "double tip-imaging," a Situation in which both tips contribute to the humelling Therefore, it has been essential to develop processes for consistently obtaining sharp, useable tips. Recently, carbon nanotubes have been used in this instance.

The hip is often mode by fungsten or plahinum

-iridium, though gold is also used. Tungsten tips are usually made by electrochemical etching, and platinum - iridium tips by mechanical shearing.

Due to the extreme sensitivity of tunnel current to height, proper vibration insulation or on extremely rigid. STM body is imperative for obtaining usable results. In the first STM by Binning and Robrer, magnetic levitation was used to keep the STM free from vibrations; now mechanical spring or gas spring systems are often used. Additionally, mechanisms for reducing eddy currents are sometimes implemented.

Maintaining the hip position with respect to the sample, scanning the sample and acquiring the data is computer controlled. The computer may also be used for enhancing the image with the help of image processing as well as performing quantitative measurements.

Controlling voltages for piezotube Piezoelectric tube with Electrodes Tunneling correct amplifier TIP Control ) AND (Scanning Sample Tunneling Voltage

'8

Operation of Scanning Tunnelling microscopy

First, a voltage bais is applied and the tip is brought close to the sample by coarse sample - to - tip control, which is turned off when the hip and sample are sufficiently dose. At dose range, fine control of the tip in all 3 dimensions when near the sample is typically piezoelectric, maintaining hip-sample seperation W typically in the 4-7 A (0.4-0.7nm) range, which is the equilibrium position between attractive (36W < 10 Å) and repulsive (W < 3Å) interactions. In this situation, the voltage bias will cause electrons to tunnel between the hip and sumple, Creating a current that can be measured. Once tunneling is established, the tip's bias and position with respect to the sample can be varied (with the details of this variation depending on the experiment) and data are obtained from the resulting changes in current.

If the hip is moved, across the sample in the x-y plane, the Changes in surface might height and density of states causes change in current. These change are mapped in images. This change in current with respect to position can be measured itself, or the height z, of the hip corresponding to a constant current can be measured. These 2 modes are called height mode and constant ament mode, respectively. In constant current mode, feedback electronics adjust the height by a voltage to the piezoelectric height control mechanism. This leads to a height variation and thus the image comes from the hip topography across the sample and gives a constant charge density Surface; this means contrast on the image is due to Variations in charge density. The benefit of using a

constant height mode is that it is foster as the plezoelectric movements require more time to register the height change in constant current mode than the current change in constant height mode. In constant height mode, the voltage and height are held both constant, while the current changes to keep the voltage from changing; this leads to an image mode of current changes over the surface, which can be related to charge density. All images produced by STM are grayscale, with colour ophionally added in past-processing in order to visually emphasize important features.

In addition to scanning across the sample, information on the electronic structure at a given location in the sample can be obtained by sweeting voltage and measuring current at a specific location. This type of measurement is called scanning tunneling spectroscopy (STS) and typically results in a plot of the local density of states as a function at energy within the sample. The advantage of stm over other measurements of the closity of states lie in it's ability to make extremely local measurement for example, the density of states at an impurity site can be compared to the density of states for from impurities.

Framerates of at least 25 Hz enable so called video-rate STM. Framerates up to 80 Hz are possible with fully working feedback that adjust the height at the hip. Due to the line-by-line scanning motion, a proter comparison on the speed requires not only the framerate, but also the number of pixels in an image, with a framerate of loke and loox 100 pixels the hip moves with a line frequency of 1 kHz, whereas it moves only with 500 Hz, when measuring with a faster framerate of Sokz but only lox 10 pixels. Video-rate STM can be used to scan surface diffusion.

guerout angen

0

aesl.in

## APPLICATIONS

OF SCANNING TUNNELING MICROSCOPY Several Surfaces have been studied with the STM. The arrangement of individual atoms on the metal surfaces on gold, platinum, nichel and copper have all been accourately documented. The absorption of and diffusion of different species such as oxygen and the epitaxial growth of gold on gold, silver on gold, and nickel on gold also have been examined in detail.

The surfaces of silicon have been studied more extensively than those of any other material. The surfaces are prepared by heated in vacuum to temperature so high that the atoms there rearrange their positions in a process called surface reconstruction. The reconstruction of the silicogn Surface designated (111) has been showed in minute detail. Such a surface constructs into an intricate and complex buttern known as the Takayang: 7x7 Structure. The position, the Chemical reactivity, and the electronic configuration of each atomic site on the 7x7 surface has been measured with the STM. The reconstruction of the silicon surface designated (100) 3 more simple. The surface atoms form pairs, or dimen, that hit into rows that extend across the entire silicon structure.

"Vacuum tunneling" of electrons from tip to sample can take place even though the environment in the region surrounding the hip is not a vacuum but is hilled with molecules of gas or liquids. With a hip-sample stowing as small as \$ 5 Å, there is little room for moleculeseven though they may exist in the surrounding atmosphere as well as in high vacuum. Indeed, it has been operated in air, in water, in insulating fluids, and In the ionic solutions used in electrochemistry. It is much more convenient than ultra-high vacuum instruments. When a high-vacuum environment is employed, its purpose is not to improve the performance of the STM but rather to ensure the cleaniness of the sample surface.

The STM can be cooled to temperatures less than 4K(-269°c or - 452°F) - the temperature of liquid Helium. It can be heated above 973K (700°C or 1300°F). The low temperature is used to investigate the properties of superconducting materials, while the high temperature is employed to shody the metals and their corrosson.

The STM is used primarily flor imaging, but there are many other modalities that have been explored. The strong electric field between the and sample has been utilized to move atoms along the sample been utilized to move atoms along the sample surface. It has been used to enhance the etching rates in various gases. In one instance, a voltage of 4 V was applied; the field at the hip was strong enough to remove atoms from the hip and deposit them on a substrate. This procedure has been employed with a gold tip to fabricate small gold structure islands or clusters on the substrate with several hundred atoms of gold in each cluster. These nanostructures are used to pattern the surface on a scale, that is unprecedented.

389950384 378

17552207